

# RFC16 SA

CATEGORY & CAPACITY  
**LAB TESTING**  
◇ CORE



**RFC16SA model is a continuous motion test tubes line to be used for preparation swab nasopharyngeal tests and testing Viruses (Ex. COVID-19).**

We can use screw-on cap (with or without pre-assembled swabs).

The machine has been designed to meet the standards in force in the pharmaceutical industry and to work in sterile rooms, plus the chance to be equipped with a LAF system.

The machine has been built using FDA-approved materials such as stainless steel and anodised aluminium. The operator can adjust most of the machine functions from the control panel.

DESCRIPTION	DATA
BOTTLE INDEX	Continuous motion
FEEDING SYSTEM	Fully automatic
MAX MECHANICAL SPEED (BPH)	8,400
WORKABLE PRODUCTS	Reagents
DOSING	Up to 10 ml of liquid product
COMPLETE AUTOMATIZATION BY TUBE HANDLING	On Optional: labeling, blistering, sleeving exc.